

Supplemental Handout for AP 2402 Lab: Hematology

Dr. Weis

Intro: Hematologic examination and testing of blood provides valuable insight regarding the health and status of the patient. This information helps the clinician in determining a diagnosis, various treatment options, and prognosis for the patient.

Purpose: The purpose of this supplemental handout is to provide students with the information needed to evaluate non-primate blood and gain experience with some of the various hematological testing methods done in a clinical lab setting. The students are directed to their text and lab manual for a more in depth discussion and presentation on blood.

Objectives :

- To safely handle all blood products and lab equipment
- To perform specific lab tests on non primate blood
- To become familiar with other testing methods for blood
- To record and analyze results from specific lab tests
- To answer discussion questions at the end of the handout

Precautions : The emphasis is SAFETY FIRST by use of gloves, cleaning of any spills and disposal of blood products and sharps in proper containers.

There is still a small chance that these animals may have contracted some unknown pathogen.

Remember :: SAFETY FIRST

GENERAL PROCEDURE :: Lab tests on Non-primate blood, outline

I. Blood smear and stain to observe

- a. RBC color, size, shape
- b. WBC differential count
- c. Platelets appearance and number

II. Hematocrit (PCV%) and Total Protein

III. Blood Typing (ABOD groups)

IV. Other tests (see lab manual for discussion)

- a. Total RBC
- b. Total WBC
- c. Hemoglobin
- d. Sedimentation rate
- e. Coagulation time

V. View other prepared slides : Leukemia, Mononucleosis, Sickle Cell Anemia

VI. Normal reference values

VII. Results Table

VIII. Questions

IX. References

Supplemental Handout for AP 2402 lab : Hematology Dr. Weis

Equipment needed :

1. Non-primate blood in proper collection tube :
 - a. Ruminant (bovine, ovine)
 - b. Equine
 - c. Canine
 - d. Feline
 - e. Rabbit
2. Slides, 1 box
3. Stain (Diff Quick / or Wrights)
4. Gloves
5. Sharps container and medical biohazard waste bag
6. Lab manual/ Text/ Handout/ Data sheet
7. Microhematocrit Centrifuge
8. Hematocrit tubes/ clay/ cardreader
9. Clinical refractometer
10. Human ABOD antisera (human antisera test kits)
11. Disposable plastic pipets
12. Toothpicks
13. Grease pen
14. Immersion oil

To use diff quick stain (modified Wright's stain), dip dried slid into each copland jar ~ 10 times at 1 second per dip. Allow excess stain to drain off before switching to next jar.

The first jar will be a light blue alcohol wash

The second jar will be a red eosin stain


The third jar will be a dark blue methylene blue stain

The fourth jar will be dH2O for the final rinse.

Please keep the lids on these stains when not in use.

Use forceps to dip the slides (or you will have stained fingers)

PCV tube
(enlarged to
show detail)



PLASMA 55%

Buffy COAT < 1%

RBC's 45%

ABO and Rh blood typing

Observation of agglutination reactions with human ABO and anti D sera in different animal species has been observed. Since the epitope for the ABO antigens consists of an oligosaccharide with four sugar components, the results can be explained due to similarities in these small saccharide chains rather than to the more variable peptide sequences of the antigens.

Results of positive agglutination reactions using animal blood and human antisera have not widely been published. From classroom observation, these are possible results * :

Canines.....B+
Rabbits.....B+
Goats.....A+
Bovine.....O+
Equine.....O+
Sheep.....O-

* remember that non-primate blood types are actually classified differently than human ABOD groups.

Mark a card or slide with the appropriate species and perform the blood typing as directed.

A positive reaction will appear as small dark spots at the margin of the mixing area. The reaction takes approximately 5 minutes. If the blood remains a homogenous pool, then you can assume that the antigen is absent.

Supplemental Handout for AP 2402 Lab : Hematology Dr. Weis
Additional Laboratory Report Questions

1. Compare and contrast mononucleosis to agranulocytic leukemia.
2. If a person with blood type A marries a person with blood type B, what are all the possible blood types for their offspring.
3. What is a transfusion reaction and why does this occur ? What are the possible consequences ?

REFERENCES

Tortora, Gerard J., Tallitsch, Robert B. Laboratory Exercises in Anatomy and Physiology, 4th ed., New York, Macmillan Publishing Company, 1993.

Evans, David L., Evans, Henriette. 1994. Using livestock blood in A&P Labs, HAPS convention handout, May 1995

Page, T.W. 1994. Blood substitutes for the physiology laboratory. HAPS News Feb., 1994

Watkins, W.M. 1966. Blood group substances in the ABO system : the genes control the arrangement of sugar residues that determine group specificity. Science 152:172.

Duncan, J. Robert, Prasse, Keith. Veterinary Laboratory Medicine Clinical Pathology, 2nd ed., Ames, Iowa State University Press, 1986.